

CE

RingCap[®]

Human Pan-Cancer Drive Gene Mutations Detection Kit

High-Throughput Sequencing

Instruction for Use



Product Name

Human Pan-Cancer Drive Gene Mutations Detection Kit (High-Throughput Sequencing)

Packing Specification

16 Tests/kit, 32 Tests/kit

Intended Use

This kit uses nucleic acids extracted from formalin-fixed paraffin-embedded (FFPE) tissues or peripheral blood plasma of cancer patients as test samples for the qualitative detection of the mutation status of relevant genes (Appendix Table 1).

The gene mutations detected by this kit cover 56 gene hotspot somatic mutations, including single-base mutations, insertions, deletions, gene fusion mutations, and other types of mutations ^[1-8]. The relevance of these mutations to targeted drugs is mainly reported from the domestic and foreign literature and has been generally accepted for clinical treatment. ^[5-10]

Technological Principles

High-Throughput Sequencing, also known as Next Generation Sequencing (NGS), can be divided into semiconductor sequencing, DNA nanosphere sequencing, and so on according to different sequencing principles. NGS enables the sequencing of up to millions of target nucleic acids at the same time, and provides abundant variation information in a short time and at a relatively low cost. Highlighting the characteristics of high output and high resolution, NGS has drawn more and more attention in multiple signaling pathways and target studies of cancer.

The construction of the sample library relies on specific modified primers and RingCap[®] mediated amplification technology with the employment of PCR apparatus. Specific modified primers enable the precise PCR amplification of target sequences, and RingCap[®] mediated amplification allows terminal modification of the products with specific sequences. With the combination of a particular PCR program and Ring-Cap[®] polymerase, library construction of target sequences could be achieved on common PCR apparatus before they are ready for high-throughput sequencing

Kit Contents

			с. ·		16 Tests/kit			32 Tests/kit		
No.	Content Name	Main Content	Color	Volume	Quantity	8-Tube Strip	Volume	Quantity	8-Tube Strip	Note
1	Onco56-DNA PCR Strip	Primer, dNTPs, Mg ²⁺ , Buffer	Blue	20 µL	16 tubes	2 tubes	20 µL	32 tubes	4 strips	Each tube contains same reagent.
2	Onco56-RNA PCR Strip	Primer, dNTPs, Mg ²⁺ , Buffer	Pink	20 µL	16 tubes	2 tubes	20 µL	32 tubes	4 strips	Each tube contains same reagent.
3	UDI 1-8 Reaction Strip	UDI primer, dNTPs, Mg ²⁺ , Buffer	Purple	20 µL	8 tubes	1 tube	20 µL	16 tubes	2 strips	Each tube represents an UDI.
4	UDI 9-16 Reaction Strip	UDI primer, dNTPs, Mg ²⁺ , Buffer	Green	20 µL	8 tubes	1 tube	20 µL	16 tubes	2 strips	Each tube represents an UDI.
5	UDI 17-24 Reaction Strip	UDI primer, dNTPs, Mg ²⁺ , Buffer	White	20 µL	8 tubes	1 tube	20 µL	16 tubes	2 strips	Each tube represents an UDI.
6	UDI 25-32 Reaction Strip	UDI primer, dNTPs, Mg ²⁺ , Buffer	Yellow	20 µL	8 tubes	1 tube	20 µL	16 tubes	2 strips	Each tube represents an UDI.
7	RingCap-Taq (1#)	Taq enzyme		20 µL	1 tube		20 µL	2 tube		
8	Onco56 Negative Control	Nuclease-Free Water		250 μL	1 tube		250 μL	1 tube		

Table1 Kit Contents



9	Onco56-DNA Positive Control	Mutation type DNA	 20 µL	1 tube	 20 µL	1 tube	
10	Onco56-RNA Positive Control	Mutation type cDNA	 20 µL	1 tube	 20 µL	1 tube	

Note 1: In UDI reaction strips, different UDI numbers respectively contain different UDI recognition sequences (see Appendix Table 3, 4). The reagents have been pre-packaged in 8-Tube strips. The left oblique position of the cap of the strip is oriented in the forward direction, from left to right followed by UDI 1, 2, 3, 4, 5, 6, 7, 8 (Figure 1).



Note 2: The contents of different batches of reagents cannot be mixed.

Additional required Equipment and Materials

- Nucleic acids extraction kit: Nucleic Acid Extraction Kit (FFPE DNA+RNA) (Xiamen Spacegen Co., Ltd, Cat. No. SPG-HSDR001R/002R) or Nucleic Acid Extraction Kit (Plasma DNA) (Xiamen Spacegen Co., Ltd, Cat. No. SPG-HSPD001R) or Nucleic Acid Extraction Kit (Peripheral blood RNA Centrifugal column method) (Xiamen Spacegen Co., Ltd, Cat. No. SPG-HSBR001R)
- 2. RNA reverses transcription kit: Super Script[™] VILO[™] cDNA Synthesis Kit (Thermo Fisher Scientific, Cat. No. 11754-050)
- 3. Quantification kit of nucleic acids: Quanti Fluor[®] dsDNA System (Promega, Cat. No. E2670) or Qubit[®] dsDNA HS Assay Kit (Thermo Fisher Scientific, Cat. No. Q32851/Q32854) or Qubit[®] ssDNA Assay Kit (Alternatively) (Thermo Fisher Scientific, Cat. No. Q10212)
- 4. Fluorometer: Quantus[™] Fluorometer (Promega, Cat. No. E6150) or Qubit[™]4.0 Fluorometer (Thermo Fisher Scientific, Cat. No. Q33238)
- 5. Magnetic beads: SG Pure Beads (Xiamen Spacegen Co., Ltd, Cat. No. SPG-PB001R/002R) or HighPrep[™] PCR (MagBio, Cat. No. AC-60005/ AC-60050/ AC-60250/ AC-60500)
- 6. Sequencing reagents and corollary reagents to be purchased separately: Selecting the corresponding sequencing reagent according to the gene sequencer
 - (1) Illumina corollary reagents: PhiX Control V3 (Illumina, Cat. No. FC-110-3001)
 - (2) MGI corollary reagents: MGIEasy universal library conversion kit (APP-A) (MGI, Cat. No. 1000004155) or High throughput sequencing primer kit (App-C) (Alternatively) (MGI, Cat. No. 1000027472)
- 7. Magnetic rack
- 8. Microvolume UV-visible spectrophotometer
- 9. Ethanol absolute (Analytical Grade)
- 10. TE Buffer (pH 8.0)
- 11. Nuclease-Free Water
- 12. Nuclease-Free pipette tips with filter

Transportation, Stability and Storage

- Storage Condition: Store the kit away from light at -15°C to -25°C, valid for 12 months. Once opened, reagents can be stored in their original packaging at -15°C to -25°C until the stated expiration date shown on the packaging. Repeated thawing and freezing should be avoided. Do not exceed a maximum of 5 freeze-thaw cycles.
- Transportation Condition: The kit should be transported at low temperature, with transporting time less than one week and transporting temperature lower than 25°C.
- 3. Check labels for the production date and expiration date of the kit.

Applicable Instruments

- 1. Library preparation PCR apparatus: ABI 9700, ABI 2720, ABI Veriti, ABI Mini Amp, etc.
- 2. Sequencing instruments:
 - (1) Illumina sequencing instruments (Miseq, NextSeq 500/550, Miniseq, etc)



(2) MGI sequencing instruments (MGISEQ-2000, DNBSEQ-G99RS, etc)

Specimen Material

The quality of the nucleic acids to be detected is critical. Please collect samples according to the following recommended sample types:

- 1. Recommended sample types: FFPE, peripheral blood plasma.
- FFPE samples: It is recommended to choose FFPE samples that have not been stored for more than 2 years and at least 30% of the collected pathological tissue is tumor lesions, and use no less than 8 pieces of 5 µm section or 5 pieces of 10 µm section for nucleic acids extraction.
- 3. Peripheral blood plasma samples: Peripheral blood should be collected with a cell-free DNA blood-collecting vessel with volumes no less than 10 mL.

Experimental Procedure

Note: Parallel library construction of Onco56-DNA Positive Control (Onco56-DNA PC), Onco56-RNA Positive Control (Onco56-RNA Positive Control (Onco56-RNA Positive Control (Onco56-RNA)) with the tested sample is suggested.

I. Library Enrichment

- 1. Reagent preparation: Thaw the **Onco56-DNA PCR Strip** (Blue) and **Onco56-RNA PCR Strip** (Pink) as needed at room temperature until no ice is present in the tubes, briefly centrifuge the strip before use. Place the **RingCap-Taq (1#)** on ice after centrifugation.
- 2. Sample preparation:
 - (1) Nucleic acid extraction and quality control:
 - (a) Commercial nucleic acids extraction kit is recommended to extract genomic DNA from the samples. Assess the quality of sample DNA with a Microvolume UV-visible spectrophotometer, the ratio of OD₂₆₀/OD₂₈₀ should be within the range of 1.8-2.2, quantify sample DNA with a Fluorometer, the concentration should be ≥2 ng/µL, the total amount of DNA should be ≥10 ng. Once the DNA quality or quantity is not conformed with the above requirements, re-extract DNA with new and/or larger input. DNA is recommended to library construction immediately or store at -15°C to -25°C for no more than 12 months.
 - (b) Commercial nucleic acids extraction kit is recommended to extract genomic RNA from the samples. Assess the quality and quantity of sample RNA with a Microvolume UV-visible spectrophotometer, the ratio of OD₂₆₀/OD₂₈₀ should be within the range of 1.8-2.3, the concentration should be ≥20 ng/µL, the total amount of RNA should be ≥100 ng. Once the RNA quality or quantity is not conformed with the above requirements, re-extract RNA with new and/or larger input. Reverse transcription is performed immediately after RNA extraction. cDNA is recommended to library construction immediately or store at -15°C to -25°C for no more than 12 months.
 - (2) DNA Sample: Dilute DNA sample to 2 ng/ μ L with TE Buffer (pH 8.0) based on the effective DNA concentration measured by the Fluorometer, and the volume $\geq 5 \mu$ L.
 - (3) cDNA Sample: cDNA sample that after reverse transcription, and the volume $\geq 5 \mu L$.
- 3. Enriching reaction for DNA
 - Add 0.25 μL of RingCap-Taq (1#) to 5 μL of the DNA Sample, Onco56-DNA PC and Onco56 NTC, vortex slightly to mix, and then centrifuge briefly.
 - (2) Gently remove the cap of the **Onco56-DNA PCR Strip**, and sequentially add 5 μ L of the template prepared above (1) into the respective tube, cap the tubes carefully.
 - (3) Centrifuge the tubes slightly and avoid creating air bubbles.
- 4. Enriching reaction for cDNA
 - Add 0.25 μL of RingCap-Taq (1#) to 5 μL of the cDNA Sample, Onco56-RNA PC and Onco56 NTC, vortex slightly to mix, and then centrifuge briefly.
 - (2) Gently remove the cap of the Onco56-RNA PCR Strip, and sequentially add 5 μL of the template prepared above (1) into the respective tube, cap the tubes carefully.
 - (3) Centrifuge the tubes slightly and avoid creating air bubbles.
- 5. Load the PCR centrifuge tubes above into the thermal cycler, then set up and run the program according to Table 2.

Table 2. PCR Amplification Procedure

Step	Temperature	Time	Cycle Number
------	-------------	------	--------------



Pre-denaturation	98°C	2 minutes	1
Denaturation	98°C	15 seconds	15
Annealing	65°C	4 minutes	15
Hold	4°C	œ	1

Note: Proceed to "Purification of Enriching Products", or store the products at 2°C to 8°C within 8 hours or at -15°C to -25°C within 24 hours. Storing for more than 24 hours is not supported.

hours. Storing for more than 24 hours is not suggested.

II. Purification of Enriching Products

Note: Transfer the magnetic beads to room temperature and vortex thoroughly to disperse magnetic beads before use. Prepare fresh 70% ethanol with Nuclease-Free Water.

- 1. Transfer 25 μ L of PCR enrichment product of **Onco56-DNA PCR Strip** each to a new 1.5 mL centrifuge tube, add 25 μ L (1×sample volume) of magnetic beads to each tube, pipet up and down 5 times to mix magnetic beads suspension thoroughly with the product.
- 2. Transfer 25 μ L of PCR enrichment product of **Onco56-RNA PCR Strip** each to a new 1.5 mL centrifuge tube, add 25 μ L (1×sample volume) of magnetic beads to each tube, pipet up and down 5 times to mix magnetic beads suspension thoroughly with the product.
- 3. Incubate the mixture for 5 minutes at room temperature.
- 4. Place the tubes on a magnetic rack for 2 minutes until the solution becomes clear, carefully remove and discard the supernatant without disturbing magnetic beads.

Note: The magnetic beads contain amplified library and should not be discarded.

- Add 150 µL of freshly prepared 70% ethanol to each tube, rotate the tubes clockwise or counterclockwise five times. Place the tubes on the magnetic rack for 2 minutes until the solution becomes clear, carefully remove and discard the supernatant without disturbing magnetic beads.
- 6. Repeat step 5 one more time for a second wash.
- 7. Remove all the ethanol from the tubes, and keep the tubes on the magnetic rack for 5 minutes at room temperature to air-dry the magnetic beads (avoid over-dry).
- 8. Remove the tube from the magnetic rack, add 35 μL of TE Buffer (pH 8.0) to each tube to fully infiltrate the magnetic beads, and vortex thoroughly (or mix by pipetting at least half the total volume up and down at least 5 times), briefly centrifuge to collect the droplets.
- 9. Incubate the mixture for 5 minutes at room temperature.
- 10. Place the tubes on the magnetic rack for 2 minutes until the solution becomes clear, carefully transfer and store the supernatant (i.e. purified product), store at -15°C to -25°C or proceed to "Library Preparation".

III. Library Preparation

Note: Use different UDI for different DNA or cDNA samples.

- 1. Reagent preparation: Thaw the **UDI Reaction Strip** based on samples amount at room temperature until no ice is present in the tubes, briefly centrifuge the strip before use. Place the **RingCap-Taq (1#)** on ice after centrifugation.
- 2. Construction reaction for DNA
 - Add 0.25 μL of RingCap-Taq (1#) to 5 μL of the purified products of DNA Sample, Onco56-DNA PC and Onco56 NTC, vortex slightly to mix, and then centrifuge briefly.
 - (2) Gently remove the cap of the UDI Reaction Strip, and sequentially add 5 μ L of the template prepared above (1) into the respective tube, cap the tubes carefully.
 - (3) Centrifuge the tubes slightly and avoid creating air bubbles.
- 3. Construction reaction for cDNA
 - Add 0.25 μL of RingCap-Taq (1#) to 5 μL of the purified products of cDNA Sample, Onco56-RNA PC and Onco56 NTC, vortex slightly to mix, and then centrifuge briefly.
 - (2) Gently remove the cap of the **UDI Reaction Strip**, sequentially add 5 μ L of the template prepared above (1) into the respective tube, cap the tubes carefully.
 - (3) Centrifuge the tubes slightly and avoid creating air bubbles.
- 4. Load the UDI Reaction Strip tubes above into the thermal cycler, then set up and run the program according to Table 3.

Table 3. PCR Amplification Procedure



Step	Temperature	Time	Cycle Number	
Pre-denaturation	98°C	2 minutes	1	
Denaturation	98°C	15 seconds	25	
Annealing	65°C	4 minutes		
Hold	4°C	x	1	

Note: Proceed to "Library Purification" or store the products at 2°C to 8°C within 8 hours or at -15°C to -25°C within 24 hours. Storing for more than 24 hours is not suggested.

IV. Library Purification

Note: Transfer the magnetic beads to room temperature and vortex thoroughly to disperse the magnetic beads before use. Prepare fresh 70% ethanol with Nuclease-Free Water.

- 1. Transfer 25 μ L of PCR product each to a new 1.5 mL centrifuge tube, add 25 μ L (1×sample volume) of magnetic beads to each tube, pipet up and down 5 times to mix the bead suspension thoroughly with the product.
- 2. Incubate the mixture for 5 minutes at room temperature.
- 3. Place the tubes on a magnetic rack for 2 minutes until the solution becomes clear, carefully remove and discard the supernatant without disturbing magnetic beads.

Note: The magnetic beads contain amplified library and should not be discarded.

- 4. Add 150 µL of freshly prepared 70% ethanol to each tube, rotate the tubes clockwise or counterclockwise five times. Place the tubes on the magnetic rack for 2 minutes until the solution becomes clear, carefully remove and discard the supernatant without disturbing magnetic beads.
- 5. Repeat step 4 one more time for a second wash.
- 6. Remove all the ethanol from the tubes, and keep the tubes on the magnetic rack for 5 minutes at room temperature to air-dry the magnetic beads (avoid over-dry).
- Remove the tubes from the magnetic rack, add 35 µL of TE Buffer (pH 8.0) to each tube to fully infiltrate the magnetic beads, and vortex thoroughly (or mix by pipetting at least half the total volume up and down at least 5 times), briefly centrifuge to collect the droplets.
- 8. Incubate the mixture for 5 minutes at room temperature.
- Place the tubes on the magnetic rack for 2 minutes until the solution becomes clear, carefully remove and store the supernatant (i.e. library) or store at -15°C to -25°C or proceed to next step.

V. Library Quantification

Bioanalyzer is recommended for the quality control of library fragments. NTC libraries should not have the fragment above 200 bp. For DNA PC library and all sample DNA libraries, the target fragments should be in 200-350 bp. For RNA PC library and all sample cDNA libraries, the target fragments should be in 200-300 bp. Fluorometer quantification kit is recommended to measure the concentration of sample library and should be more than 0.5 ng/µL.

VI. Library dilution, mixing and sequencing

- 1. Illumina sequencing platform
 - (1) According to the library concentration measured by the Fluorometer, use the following formula to convert the molar concentration of the library.

Library concentration (nM) = $\frac{\text{Library concentration } (ng/\mu L) \times 10^{6}}{\text{Library length } (bp) \times 650}$

- (2) Per the concentration measured, dilute the sample library to 4 nM with Nuclease-Free Water.
- (3) The proportion of DNA and cDNA is 20:1 (Mix 20 µL of each DNA sample library with 1 µL of cDNA sample library).
- (4) Proceed sample dilution and denaturation according to the matching Illumina sequencing kit (refer to the operation manual of each equipment).
- (5) The concentration of Phix Control V3 is more than 5% (for example: If the loading volume is 600 μ L, the volume occupied by Phix Control V3 should be more than 30 μ L).
- (6) Library on-machine sequencing (according to the specification of instrument and matching reagent).
- 2. MGI sequencing platform
 - (1) The recommended amount of cyclic library input is 0.5 pmol, the loading ratio of DNA library and cDNA library is 20:1, and the



required proportion of each library in the total library of 0.5 pmol is calculated, and then the required amount of each library is calculated according to the following formula:

Library input (ng)= Library length (bp) × 650×proportion%

1000

- (2) According to the measured library concentration, the required input volume is calculated, and then mixed to obtain 0.5 pmol total library, the total volume is not more than 48 μL.
- (3) Denaturation and cyclization of libraries according to the MGIEasy Universal Library Conversion Kit (App-A) (no need for splitconversion PCR, see the accompanying reagent manual for instructions).
- (4) Proceed sample dilution and denaturation according to the matching sequencing kit (refer to the operation manual of each equipment).
- (5) Library DNB preparation and on-machine sequencing (operation instructions refer to the accompanying reagent manual and instrument manual).

Note: Store undiluted libraries at -15°C to -25°C for up to 7 days, the mixture of diluted libraries is suggested to be used right after it is ready.

VII. Bioinformatics Analysis

Transfer the Fastq files obtained by sequencing to the analysis server, perform data quality control, sequence alignment, mutation annotation and gene fusion analysis-based on the Clinical NGS Data Analysis System of Xiamen Spacegen Co., Ltd.

Data Analysis

- 1. Results of DNA
 - Standard of quality: For all sample DNA libraries, the target fragment should be in 200-350 bp, On Target should be ≥80%, Uniformity should be ≥75% and mean depth should be ≥5000×.
 - (2) Mutated positive judgement criterion: In the result of variations analysis, if effective depth is >500×. And mutation frequency is >1%, this mutation site is judged as positive mutation. Otherwise, it is judged as negative or below the detection limit.
- 2. Results of RNA
 - (1) Standard of quality: For all sample cDNA libraries, the target fragments should be in 200-300 bp. The value of "Total Reads" of at least 2 of the 5 internal control genes (HMBS, TBP, MYC, LRP1, and MRPL13) recorded in the file should be ≥200×, that guarantees the quality of RNA sample.
 - (2) Mutated positive judgement criterion: The interpretation mode of fusion mutation provided by this kit is "specific fusion site detection". The sequence data obtained from fusion transcription analysis is recorded in the format of several readings per target, which is different from other analysis platforms. The background signal level should be taken into account when interpreting the results:
 - (a) If the sample is not read in both positive and negative directions, indicate the specific fusion site is negative.
 - (b) If the sample is read in both positive and negative directions, and the Total Reads is less than 200×, indicates that the fusion close to the background signal has been detected, it is recommended to increase the input of the library; After re-detection, the Total Reads is still less than 200× indicate the specific fusion site is negative or lower than the detection limit.
 - (c) If the sample is read in both positive and negative directions, and Total Reads ≥200×, indicate the specific fusion site is positive.

Interpretation of Results

- 1. For the negative control library may be ≥ 0 ng/ μ L, but it should not have any fragment above 200 bp. Otherwise, this test is invalidated.
- For the DNA positive control library, the target fragment should be in 200-350 bp, as well as On Target should be ≥80%, Uniformity should be ≥75%, moreover, mean depth ≥5000×. Otherwise, this test is invalidated.
- 3. For the RNA positive control library, the target fragments should be in 200-300 bp, it should more than 2 of the 5 internal control genes are all read, and the Total Reads \geq 200×. Otherwise, this gene fusion test results are invalidated.
- 4. For the DNA sample libraries, the target fragment should be in 200-350 bp, each amplicon should have coverage, as well as On Target should be ≥80%, Uniformity should be ≥75%, moreover, mean depth should be ≥5000×. Otherwise, this mutation detection results are invalidated.
- 5. For the RNA sample libraries, the target fragments should be in 200-300 bp, it should more than 2 of the 5 internal control genes are all read, and the Total Reads \geq 200×. Otherwise, this gene fusion test results are invalidated.
- 6. The grade of somatic variation based on the "Standards and Guidelines for the Interpretation and Reporting of Sequence Variants in



Cancer: A Joint Consensus Recommendation of the Association for Molecular Pathology, American Society of Clinical Oncology, and College of American Pathologists" jointly formulated by AMP/ASCO/CAP in 2017 could divide into 4 types:

- Clear clinical significance: Diagnostic/prognostic marker of specific tumor or drugs recommended/approved in the professional guidelines.
- (2) Potential clinical significance: Diagnostic prognostic marker of specific tumor or drugs that has level A evidence of another tumor in the multiple small research.
- (3) Unknown clinical significance: It is not found higher rates of variants in the general population and tumor databases, moreover, not has clear published evidence.
- (4) Harmless or may be harmless clinical significance: It is found higher rates of variants in the general population and not published evidence.

Limitation of the Kit

- 1. The detection results are for research use only. For mutation sites that are not included in the kit, or the nucleic acids extracted from samples are stored longer than required, the results shall not be interpreted by the instruction.
- 2. The negative results cannot exclude the mutations. For few tumor cells, excessive degradation, or the nucleic acids concentration is below the detection limit can also cause a negative result.
- 3. Unreasonable sample collection, transportation, processing, improper operation and the experimental environment may lead to false negative or positive results.
- 4. Tumor tissue (cells) may have large heterogeneity, different test results may be obtained by sampling different parts.

Performance characteristics

- 1. The kit should be neat in appearance, clearly labels, and no leakage.
- 2. When unfrozen, the reagents shall be clear, without sediments.
- 3. Negative reference conformity rate should be 100%.
- 4. Positive reference conformity rate should be 100%.
- 5. The kit allows the detection of 1% of specific gene mutations in 10 ng DNA tissue samples.
- 6. The kit allows the detection of 20 copies/ μ L of fusion mutations in RNA tissue samples.

Warnings and Precautions

- 1. Please read the instruction carefully in prior to experiments.
- 2. Conduct experiments abided by laboratory regulations to reduce cross-contaminations of products or reagents; divide experiment areas into different function zones if possible.
- 3. Avoid repetitively freezing and thawing the reagents in the kit. Do not exceed a maximum of 5 freeze-thaw cycles.
- 4. The results of this kit will be affected by sample source, collection process, quality, transportation conditions, pre-treatment, etc., as well as the quality of the extracted nucleic, instrument types, operating environment, and the limitation of current molecular biotechnology. The factors and variables mentioned above would lead to false positive or false negative results. Users must be aware of the potential errors, accuracy and limitations that may exist during the process of detection.
- 5. The quality of nucleic acids is crucial, and the quality control of DNA should be performed after extraction, proceed to further steps immediately or store properly at -15°C to -25°C. RNA is recommended to be reverse transcript to cDNA before storage, and RNA without reverse transcription is recommended to be stored below-70°C.
- 6. Do not substitute any original reagents contained in the kit. Do not mix reagents with different Lots.
- 7. Pay special attention to the use of positive control and the use of filter pipette tips is highly recommended to avoid false-positive results caused by contamination of reagents.
- 8. Be cautious of contamination from external nucleic. Ensure to add the nucleic template before operating the positive control. Segregate areas for reagent preparation and sample processing. Use dedicated pipettes and pipette tips for reagent preparation and template addition, respectively.
- 9. Clean experiment areas before experiment with 10% hypochlorous acid followed by twice water rinsing. Sterilize the environment and pipettes with 10% hypochlorous acid, 75% ethanol, or UV radiation after experiment.
- 10. All reagents in use have potential hazard. It is recommended wearing proper protective suit and gloves. For first-use of this kit, you may



receive training by our technical supports.

- 11. All samples including positive control in the kit should be considered as potential infectious substances which should be handled carefully.
- 12. Avoid using peripheral wells of PCR instrument; vacate holes or columns between samples to avoid cross-contamination.

Symbols

Symbol	Symbol definition
	Indicates the need for the user to consult the instructions for use.
IVD	Indicates a medical device that is intended to be used as an in vitro diagnostic medical device.
	Indicates the date when the medical device is manufactured.
LOT	Indicates the manufacturer's batch code so that the batch or lot can be identified.
X	Indicates the temperature limitation.
	Indicates the date after which the medical device is not to be used.
$\underbrace{\uparrow} \underbrace{\uparrow}$	This is the correct upright position of the distribution packages for transport or storage.
Ţ	Indicates a medical device should be kept dry.
溇	Indicates a medical device that needs protection from light sources.
	Indicates the medical device manufacturer.
EC REP	Indicates the authorized representative in the European Community/European Union.
()	The product meets the basic requirements of European in vitro diagnostic medical devices directive 98/79/EC.

References

- MorkTs, Wu YL, Thongprasert S, et al. Gefitinib or carboplatin-paclitaxel in pµlmonaryademocarcinoma. N Engl J Med.2009,361 (10):947-57.
- 2. Gazdar AF. Personalized medicine and inhibition of EGFR singnaling in lung cancer. N Engl J Med. 2009, 361 (10): 1018-20.
- 3. Dancey JE. Epidermal growth factor receptor inhibitors in non-small cell lung cancer. Drugs. 2007, 67 (8): 1125-38.
- 4. Jemal A, Bray F, Center VI M, el al. Global cancer statistics.CA Cancer J Clin. 2011. 61 (2): 69-90.
- 5. Souliares D,Greer W,Maglioeco AM,et al. KRAS mutationtesting in the treatment of metastatic colorectal cancer with anti-EGFR



therapies. CurrOncol.2010,17 Suppl 1: S31-S40.

- 6. Paul B.Chapman, M.D., et al. Improved Survival with Vemurafenibin Melanoma with BRAF Mutation. N Engl J Med. 2011:364,2507-16.
- 7. Wendy De Roock, Bart Claes, et al. Effects of KRAS, BRAF, NRAS, and PIK3CA mutationson the efficacy of cetuximab plus chemotherapy inchemotherapy-refractory metastatic colorectal cancer: a retrospective consortium analysis. Lancet Oncol. 2010, 11:753–62.
- 8. Soda M,Choi YL, Enomoto M,et al.Identification of the transforming EML4-ALK fusion gene in non-small-cell lung cancer. Nature.2007, 448 (7153): 561-6.
- 9. Sai-Hong Ignatius Ou, Jackie Tan, et al. ROS1 as a 'druggable' receptor tyrosine kinase: lessons learned from inhibiting the ALK pathway. Expert Rev. Anticancer Ther.2012,12 (4):447-56.
- Alexa B.Turke, Kreshnik Zejnullahu, et al. Preexistence and Clonal Selection of MET Amplification in EGFR Mutant NSCLC. Cancer Cell. 2010, 17: 77-88.

Г 7



Lotus NL B.V.

Address:Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands. E-mail: peter@lotusnl.com

Manufacturer: XIAMEN SPACEGEN CO., LTD. Address: 4th floor, No.2041 Xizhou Road, Tong'an District, Xiamen 361100, P. R. China Tel: +86 592 7578317 Fax: +86 592 7578319 E-mail: spacegen@ispacegen.com Website: http://www.sspacegen.com



Appendix Table 1. Gene Information Included in the Kit									
DNA-Gene									
ABL1	AKT1	ALK	APC	ATM	BRAF	CDH1	CDKN2A		
CSF1R	CTNNB1	EGFR	ERBB2	ERBB4	EZH2	FBXW7	FGFR1		
FGFR2	FGFR3	FLT3	GNA11	GNAQ	GNAS	HNF1A	HRAS		
IDH1	IDH2	JAK2	JAK3	KDR	KIT	KRAS	MET		
MLH1	MPL	NOTCH1	NPM1	NRAS	PDGFRA	PIK3CA	PTEN		
PTPN11	RB1	RET	SMAD4	SMARCB1	SMO	SRC	STK11		
TP53	VHL								

 RNA(Fusion)-Gene

 ALK
 ROS1
 RET
 NTRK1
 NTRK3
 MET exon14- skipping

Positive Control	Gene	Base Mutation	Amino Acid Mutation	Cosmic ID	Mutation Type
	EGFR	c.2235_2249del15	p.E746_A750delELREA	6223	Deletion mutation
Onco56 DNA	KRAS	c.35G > A	p.G12D	521	Point mutation
Positive Control	BRAF	c.1799T > A	p.V600E	476	Point mutation
	HER2 c		p.A775_G776insYVMA	20959	Insertion mutation
Onco56 RNA	ALK	EML4-ALK(E13-A20)		463	Fusion mutation
Positive Control	ROS1	SLC34A2-ROS1(S4-R32)		1197	Fusion mutation

Appendix Table 3. Information of UDI Recognition Sequences based on Illumina Tech

	NT 1	.7.0	i5 Sequence	i5 Sequence	
Strip Color	Number	1/ Sequence	(NovaSeq, MiSeq)	(iSeq, MiniSeq, NextSeq)	
	UDI-1	TGCATAGC	TAGGATTC	GAATCCTA	
	UDI-2	TCTATGCA	GTCGTTGC	GCAACGAC	
	UDI-3	GTACGCAT	CCTCGCAT	ATGCGAGG	
D	UDI-4	AGGTCCTG	AGAAGGCG	CGCCTTCT	
Purple	UDI-5	CATGAGCT	ACGTCAGA	TCTGACGT	
	UDI-6	AACTCTAG	CATCTGAT	ATCAGATG	
	UDI-7	CCGGATGC	GTATCACG	CGTGATAC	
	UDI-8	GTACGATA	TGCAACTA	TAGTTGCA	
	UDI-9	ATTCGATA	ATGGATCG	CGATCCAT	
	UDI-10	CGTAGTAC	GCTGAATG	CATTCAGC	
	UDI-11	GAGTACGT	CAACTGGC	GCCAGTTG	
Cassa	UDI-12	TCAGTGCG	TGCAGCAT	ATGCTGCA	
Green	UDI-13	CACACAGT	ACGACCAA	TTGGTCGT	
	UDI-14	GTGCATCG	CATTCGGC	GCCGAATG	
	UDI-15	TGCGTCAC	GTATGATT	AATCATAC	
	UDI-16	ACATCGTA	TGCCTTCA	TGAAGGCA	
	UDI-17	CGGAACGA	GCTGGCTT	AAGCCAGC	
	UDI-18	CCTGGCAC	ATAGAGAC	GTCTCTAT	
	UDI-19	ATATCGCT	CACATTGA	TCAATGTG	
White	UDI-20	GACAGTTG	TGGTCACG	CGTGACCA	
white	UDI-21	TGCCTATG	ACCTTCGG	CCGAAGGT	
	UDI-22	GTACCAGT	CGACCATC	GATGGTCG	
	UDI-23	AATGTGCA	TAGCATCA	TGATGCTA	
	UDI-24	TCGTATAC	GTTAGGAT	ATCCTAAC	
	UDI-25	CTGTGTGT	CGTCGTCT	AGACGACG	
	UDI-26	ACAGCACT	ATCCTAGC	GCTAGGAT	
	UDI-27	TATCAGTG	GAAGCCTG	CAGGCTTC	
Vallow	UDI-28	CGGTGTTA	TCGAAGTA	TACTTCGA	
Tenow	UDI-29	GTCATCAC	ACCGGTAC	GTACCGGT	
	UDI-30	GATGTCAG	CATTCAAT	ATTGAATG	
	UDI-31	TCACAGCA	TGGTAGCA	TGCTACCA	
	UDI-32	AGCACAGC	GTAATCGG	CCGATTAC	

Appendix Table 4. Information of UDI Recognition Sequences based on MGI Tech

Strip Color	Number	Sequence		Strip Color	Number	Sequence
	UDI-1	TAGGATTCTGCATAGC GTCGTTGCTCTATGCA CCTCGCATGTACGCAT AGAAGGCGAGGTCCTG		White	UDI-17	GCTGGCTTCGGAACGA
	UDI-2				UDI-18	ATAGAGACCCTGGCAC
Purple	UDI-3				UDI-19	CACATTGAATATCGCT
	UDI-4				UDI-20	TGGTCACGGACAGTTG
	UDI-5	ACGTCAGACATGAGCT			UDI-21	ACCTTCGGTGCCTATG



	1	1	1	1	1	
	UDI-6	CATCTGATAACTCTAG			UDI-22	CGACCATCGTACCAGT
	UDI-7	GTATCACGCCGGATGC			UDI-23	TAGCATCAAATGTGCA
	UDI-8	TGCAACTAGTACGATA			UDI-24	GTTAGGATTCGTATAC
	UDI-9	ATGGATCGATTCGATA			UDI-25	CGTCGTCTCTGTGTGT
	UDI-10	GCTGAATGCGTAGTAC			UDI-26	ATCCTAGCACAGCACT
	UDI-11	CAACTGGCGAGTACGT			UDI-27	GAAGCCTGTATCAGTG
Graan	UDI-12	TGCAGCATTCAGTGCG		Yellow	UDI-28	TCGAAGTACGGTGTTA
Green	UDI-13	ACGACCAACACACAGT			UDI-29	ACCGGTACGTCATCAC
	UDI-14	CATTCGGCGTGCATCG			UDI-30	CATTCAATGATGTCAG
	UDI-15	GTATGATTTGCGTCAC			UDI-31	TGGTAGCATCACAGCA
	UDI-16	TGCCTTCAACATCGTA			UDI-32	GTAATCGGAGCACAGC