



RingCap[®]

BRCA1/2 Gene Mutations Detection Kit

High-Throughput Sequencing

Instruction for Use (Illumina)

Product Name

BRCA1/2 Gene Mutations Detection Kit (High-Throughput Sequencing)

Packing Specification

16 Tests/Kit, 32 Tests/Kit

Intended Use

The kit is intended for the detection of whole exon and exon-intron junctions in BRCA1/2 with FFPE pathological tissue and peripheral blood samples from patients with breast cancer or ovarian cancer. The results are only for scientific reference.

BRCA1/2 is the abbreviation of Breast Cancer Susceptibility Gene 1/2; Both of them are tumor suppressor genes and play important roles in DNA damage repair, cell cycle regulation, gene transcription activation, and chromatin stability. BRCA1 gene(NM_007294) locates in human chromosome 17q21, is composed of 23 exons; BRCA2 gene(NM_000059) locates in human chromosome 13q12.3, is composed of 27 exons.

Among women suffer from breast cancer and ovarian cancer, about 2% and 10-15% respectively of the patients have been reported to carry the mutations of BRCA1/2, moreover, the mutation status is closely related to the familial breast cancer, namely, people with BRCA1/2 mutations may be more susceptible to breast cancer and/or ovarian cancer. It is studied that people who're mutation carriers may suffer from breast cancer at the probability of 40-80% and ovarian cancer at 16-60%, while the probability is 12% and 1% respectively among people with wild-type phenotype. Furthermore, BRCA2 mutation may also lead to the occurrence of breast cancer in male with the occurrence rate at about 6%.

Neither BRCA1 nor BRCA2 gene has mutations of high frequency, which means that hotspot sites sequencing is inadequate, but call for complete genome sequencing then. This kit facilitates library construction of whole exome of BRCA1/2 genes before quick detection of 2229 mutation sites (see Attached Table 1) with the assist of high throughput sequencer. Scientific and clinical trials have shown that the new drug, olaparib, could be a target therapy for cancer patients with BRCA1 or BRCA2 mutation. The detection of BRCA1/2 gene mutation status has gradually become an aid of individualized therapy in clinic.

Technological Principle

High-Throughput Sequencing, also known as Next Generation Sequencing (NGS), can be divided into semiconductor sequencing, DNA nanosphere sequencing and so on according to the different sequencing principles. NGS enables the sequencing of up to millions of target nucleic acids at once, provides abundant variation information in a short time and at a relatively low cost. Highlighting the characteristics of high output and high rebuffer, NGS has drawn more and more attention in multiple signaling pathways and targets studies of cancer.

The construction of sample library relies on specific modified primers and RingCap® mediated amplification technology with the employment of PCR apparatus. Specific modified primers enable the precise PCR amplification of target sequences, and RingCap® mediated amplification allows terminal modification of the products with specific sequences. With the combination of a particular PCR program and Ring-Cap® enzyme, library construction of target sequences could be achieved on common PCR apparatus before they are ready for high-throughput sequencing.

Kit Contents

Table 1. Kit Contents

No.	Content Name	Main Content	Strip Color	16 Tests/Kit			32 Tests/Kit			Note
				Volume	Quantity	8-Tube Strip	Volume	Quantity	8-Tube Strip	
1	ILL-BRCA enriching PCR strip 1	Primer, dNTPs, Mg ²⁺ , buffer	Blue	20 μL	16 tubes	2 strips	20 μL	32 tubes	4 strips	Each tube contains same reagent.

2	ILL-BRCA enriching PCR strip 2	Primer, dNTPs, Mg ²⁺ , buffer	Pink	20 μL	16 tubes	2 strips	20 μL	32 tubes	4 strips	Each tube contains same reagent.
3	Index 1-8 ligation reaction strip	Barcode, dNTPs, Mg ²⁺ , buffer	Purple	20 μL	8 tubes	1 strip	20 μL	8 tubes	1 strip	Each tube represents a barcode.
4	Index 9-16 ligation reaction strip	Barcode, dNTPs, Mg ²⁺ , buffer	Green	20 μL	8 tubes	1 strip	20 μL	8 tubes	1 strip	Each tube represents a barcode.
5	Index 17-24 ligation reaction strip	Barcode, dNTPs, Mg ²⁺ , buffer	White	—	—	—	20 μL	8 tubes	1 strip	Each tube represents a barcode.
6	Index 25-32 ligation reaction strip	Barcode, dNTPs, Mg ²⁺ , buffer	Yellow	—	—	—	20 μL	8 tubes	1 strip	Each tube represents a barcode.
7	RingCap-Taq (1#)	Taq enzyme	—	20 μL	1 tube	—	20 μL	2 tubes	—	—
8	BRCA Negative Control	Purified water	—	1.0 mL	1 tube	—	1.0 mL	1 tube	—	—
9	BRCA Positive Control	Wild-type cell line DNA	—	20 μL	1 tube	—	40 μL	1 tube	—	—

Note: In Index reaction strips, different index numbers respectively contain 32 different IIIIdx recognition sequences (see Appendix Table 1); the reaction buffer has been pre-loaded in the 8-tube strips; the left oblique position of the cap of the strip is oriented in the forward direction, from left to right followed by index 1, 2, 3, 4, 5, 6, 7, 8 (Figure 1).




Figure 1. Index numbers of 8-tube strips

Note: The components of different batches of reagents cannot be mixed.

Equipment and Reagents Required

1. Fluorometer: Promega Quantus™ Fluorometer, Cat. No. E6150; Thermo Fisher Scientific Qubit® 4.0 Fluorescence Meter, Cat. No. Q32866;
2. Nucleic acid extraction kit: commercial nucleic acid extraction kits are recommended;
3. Quantification kit of nucleic acids: Promega, QuantiFluor® dsDNA System, Cat. No. E2670; Thermo Fisher Scientific Qubit® dsDNA HS Assay Kit, Cat. No. Q32851/Q32854;
4. Magnetic beads: Xiamen Spacegen Co., Ltd SGpure beads, Cat. No. SPG-PB001; Magnetic beads Kit from Beckman Coulter, Cat. No. A63880/A63881/A63882;
5. Sequencing reagents: Selecting the corresponding sequencing reagent according to the gene sequencer;
6. TE buffer (pH 8.0);
7. Absolute ethanol (Analytical Grade);
8. Nuclease-free water;
9. Nuclease-free pipettes and tips.

Storage and Stability

1. Storage Condition. Store the kit away from light at -20±5°C, valid for 9 months, and is not influenced by bottle openings, moreover, do not use the reagents after 5 repeated freeze-thaw cycles.
2. Transportation Condition. The kit should be transported in foam cases with ice bags, with transporting time of less than one week and transporting temperature lowers than 25 °C.
3. Check the labels for the production date and expiration date of the kit.

Applicable Instruments

1. PCR system/ thermal cycler: ABI9700, ABI 2720, ABI Veriti, ABI MiniAmp;
2. Sequencing instruments: Miseq, NextSeq 500/550, Miniseq.

Specimen Material

The quality of the DNA is critical. Therefore, collecting samples according to the following recommended sample types and requirements, followed by DNA extraction:

1. Recommended sample types: FFPE pathological tissue or slices, peripheral blood;
2. FFPE samples: ensure that at least 20% of the collected pathological tissue is tumor lesions; choose FFPE samples that have not been stored for more than 2 years; extract DNA with at least 8 slices of 5 μm section or at least 5 slices of 10 μm section;
3. Peripheral blood: Peripheral blood should be collected with a non-heparin anticoagulant with volumes no less than 2 mL;
4. Commercial kits are highly recommended to extract genomic DNA from the samples. Quantify sample DNA with a fluorescence meter, the concentration should be $\geq 2 \text{ ng}/\mu\text{L}$, total amount of DNA should be $\geq 20 \text{ ng}$. Once the DNA quantity or quality did not conform to the above requirements, re-extract DNA with a resample or a larger amount of samples. Proceed to library construction after DNA extraction or store at $-20\pm 5 \text{ }^\circ\text{C}$ for less than 12 months.

Experimental Procedure

Note: Parallel library construction of BRCA Positive Control (PC) and BRCA Negative Control (NTC) with tested sample is suggested.

I. Library Enrichment

1. Reagent preparation: unfreeze the **ILL-BRCA enriching PCR strip 1 (blue)** and **ILL-BRCA enriching PCR strip 2 (pink)** at room temperature according to the total number of samples, briefly centrifuge the tubes before use; place the **RingCap-Taq (1#)** on ice after centrifugation;
2. Sample preparation: dilute sample DNA to 2 $\text{ng}/\mu\text{L}$ with TE buffer (pH 8.0), and prepare $\geq 10 \mu\text{L}$ of the diluted sample;
3. Enriching PCR reaction:
 - a) Pipet 0.5 μL RingCap-Taq (1#) to 10 μL of the “DNA Sample”, “BRCA-PC”, and “BRCA-NTC”, vortex slightly followed by brief centrifugation;
 - b) Gently remove the cap of enriching PCR strip, for one sample, respectively pipet 5 μL of the template prepared above (10 μL) into ILL-BRCA enriching PCR tube 1 (blue) and ILL-BRCA enriching PCR tube 2 (pink), replace the cap carefully;
 - c) Centrifuge the tubes/ strips slightly to dislodge bubbles;
4. Load the PCR reaction tubes / strips into the PCR instrument, Set up and perform the amplification procedure according to Table 2:

Table 2. PCR amplification procedure

Step	Temperature	Time	Cyclic number
Pre-denaturation	98 $^\circ\text{C}$	2 minutes	1
Denaturation	98 $^\circ\text{C}$	15 seconds	15
Annealing	65 $^\circ\text{C}$	4 minutes	
Storage	4 $^\circ\text{C}$	∞	1

Note: Proceed to “Purification of Enriching Products”, or store the products at 2 ~ 8 $^\circ\text{C}$ within 8 hours or at $-20\pm 5 \text{ }^\circ\text{C}$ within 24 hours.

Storing for more than 24 hours is not suggested.

II. Purification of Enriching Products

Note: Bring magnetic beads to room temperature and vortex thoroughly to disperse magnetic beads before use; prepare fresh 70% ethanol with nuclease-free water.

1. Mix the two PCR enrichment products of one sample into a new 1.5 mL eppendorf tube, add 50 μL magnetic beads to each tube, pipet up and down to mix the bead suspension thoroughly with the product;
2. Incubate the mixture for 5 minutes at room temperature;
3. Place the tube on a magnet rack, then incubate for 2 minutes, carefully remove and discard the supernatant without disturbing magnetic beads;
4. Pipet 150 μL of freshly prepared 70% ethanol into each tube, rotate the tube side-to-side in the two positions of the magnetic rack 5 times to wash magnetic beads, place the tube on the magnetic rack for 2 minutes, carefully remove and discard the supernatant without

- disturbing magnetic beads;
- Repeat step 4 for a second wash;
 - Remove all the ethanol from the tube, and keep the tube on a magnetic rack for 5 minutes to air-dry magnetic beads (avoid over-dry);
 - Remove the tube from the magnetic rack, add 35 μL of TE buffer (pH 8.0) to each tube, and vortex thoroughly (alternatively, mix by pipetting at least half the total volume up and down), briefly centrifuge to collect the droplets, Incubate the mixture for 5 minutes at room temperature;
 - Place the tube on a magnetic rack for 2 minutes until the solution is clear, carefully remove and store the supernatant (i.e. **purified product**) at $-20\pm 5^\circ\text{C}$ or proceed to “Library Construction”.

III. Library Construction

- Reagent preparation: unfreeze the **Index ligation reaction strip** at room temperature, briefly centrifuge the tubes before use; place the **RingCap-Taq (1#)** on ice after centrifugation;
- BRCA library construction:**
 - Pipet 0.25 μL **RingCap-Taq (1#)** to 5 μL of the purified products of “DNA Sample”, “BRCA-PC”, and “BRCA-NTC”, vortex slightly followed by brief centrifugation;
 - Gently remove the cap of the **index ligation reaction strip**, sequentially pipet 5 μL of the template prepared above into the respective tube, and replace the cap carefully;
 - Centrifuge the tubes slightly to dislodge bubbles;
- Load the PCR reaction tubes/ strips into the PCR instrument, and set up and perform the amplification procedure according to Table 3:

Table 3. PCR amplification procedure

Step	Temperature	Time	Cyclic number
Predenaturation	98 $^\circ\text{C}$	2 minutes	1
Denaturation	98 $^\circ\text{C}$	15 seconds	25
Annealing	65 $^\circ\text{C}$	4 minutes	
Storage	4 $^\circ\text{C}$	∞	1

Note: Proceed to “Library Purification”, or store the products at $2\sim 8^\circ\text{C}$ within 8 hours or at $-20\pm 5^\circ\text{C}$ within 24 hours. Storage for more than 24 hours is not suggested.

IV. Library Purification

Note: Bring magnetic beads to room temperature and vortex thoroughly to disperse magnetic beads before use; prepare fresh 70% ethanol with nuclease-free water.

- Transfer 25 μL of the PCR library product to a new 1.5 mL Eppendorf tube, add 25 μL of magnetic beads to each tube, and pipet up and down to mix the bead suspension thoroughly with the product;
- Incubate the mixture for 5 minutes at room temperature;
- Place the tube on a magnet rack, then incubate for 2 minutes, carefully remove and discard the supernatant without disturbing magnetic beads;
- Pipet 150 μL of freshly prepared 70% ethanol into each tube, rotate the tube side-to-side in the two positions of the magnetic rack for 5 times to wash magnetic beads, place the tube on a magnetic rack for 2 minutes, carefully remove and discard the supernatant without disturbing magnetic beads;
- Repeat step 4 for a second wash;
- Remove all the ethanol from the tube, and keep the tube on a magnetic rack for 5 minutes to air-dry magnetic beads (avoid over-dry);
- Remove the tube from magnetic rack, add 35 μL of TE buffer (pH 8.0) to each tube, and vortex thoroughly (alternatively, mix by pipetting at least half the total volume up and down), briefly centrifuge to collect the droplets; Incubate the mixture for 5 minutes at room temperature;
- Place the tube on a magnetic rack for 2 minutes until the solution is clear, carefully remove and store the supernatant (i.e. **library**) at -

20±5°C or proceed to “Library Quantification and Dilution”.

V. Library Quantification and Dilution

1. Quality control (QC) of sample library: A bioanalyzer is recommended for the quality control of library fragments; for the BRCA-NTC library, no fragments shall be detected; for BRCA-PC and all sample libraries, the main fragments should be at 250 ~ 350 bp; for library effective concentration lowers than 1 ng/μL (measured by fluorometer), either of which is decided unqualified;
2. The concentration of Phix Control V3 more than 5% (e.g.the percentage of Phix Control V3 should be more than 30 μL in the 600 μL loading volume);
3. Sample dilution and denaturation according to the matching Illumina sequencing kit (refer to the operation manual of each piece of equipment);
4. Store undiluted sample libraries at - 20±5°C for up to 7 days; The mixture of diluted libraries is suggested to be used right after it is ready.

VI. Bioinformatics Analysis

Transfer the Fastq files obtained by sequencing to the analysis server, followed perform data quality control, sequence alignment, mutation annotation analysis-based on the clinical high-throughput sequencing data analysis system (abbreviated as analysis system below) of Xiamen Spacegen Biotechnology Co., Ltd.

Positive Judgment Value

1. Standard quality control: For DNA sample libraries, the main fragments should be at 250 ~ 350 bp, Uniformity ≥ 75%, mean Depth ≥ 500×;
2. Mutated positive judging criteria:
 - a) Peripheral blood sample: In the result of germline variation analysis, if effective depth > 100× and mutation frequency > 25%, this mutation site is judged as positive mutation; Otherwise, it is judged as negative.
 - b) FFPE: In the result of somatic variation analysis, if effective depth > 100× and mutation frequency > 5%, the site of this mutation is judged as positive mutation; Otherwise, it is judged as negative or below the detection limit.

Interpretation of Results

1. NTC libraries should not any fragment; Otherwise, this test is invalidated;
2. For DNA positive control library, the target fragment should be in 250 ~ 350 bp as well as Uniformity should be more than 75%, moreover, Mean Depth more than 500×;
3. For DNA tested sample library, the target fragment should be in 250 ~ 350 bp as well as Uniformity should be more than 75%, moreover, Mean Depth more than 500×; Otherwise, this test is invalidated;
4. The grade of somatic variation based on the "Standards and Guidelines for the Interpretation and Reporting of Sequence Variants in Cancer: A Joint Consensus Recommendation of the Association for Molecular Pathology, American Society of Clinical Oncology, and College of American Pathologists" jointly formulated by AMP/ASCO/CAP in 2017 could divide into 4 types:
 - a) Clear clinical significance: Diagnostic\prognostic marker of specific tumor or drugs recommended\approved in the professional guidelines;
 - b) Potential clinical significance: Diagnostic\prognostic marker of specific tumor or drugs that has level A evidence of another tumor in the multiple small research;
 - c) Unknown clinical significance: It is not found higher rates of variants in the general population and tumor databases, moreover, not has clear published evidence;
 - d) Harmless or may be harmless clinical significance: It is found higher rates of variants in the general population and not published evidence.
5. The grade of germline variation based on the "Standards and guidelines for the interpretation of sequence variants: a joint consensus recommendation of the American College of Medical Genetics and Genomics and the Association for Molecular Pathology" and another standard of genetic variation, variation types, and conserved biological functional prediction, evidence of databases(e.g.ClinVar、COSMIC) and literature could divide into 5 levels
 - a) Pathogenic: There is clear evidence that this variant is associated with an increased risk of cancer;
 - b) Likely pathogenic: There is moderate evidence that this variant is associated with an increased risk of cancer;

- c) Uncertain significance: There is unclear evidence that this variant is associated with an increased risk of cancer;
- d) Likely benign: There is moderate evidence that this variant is unrelated to an increased risk of cancer;
- e) Benign: There is clear evidence that this variant is unrelated to an increased risk of cancer.

Limitation of the Kit

Testing results obtained from the kit should only be taken as a scientific reference. The instruction shall not interpret the results for mutation sites that were not included in the kit or DNA extracted from samples that were not collected according to designated requirements.




Physical Performance of Products










1. The kit should have a neat appearance, clear markings, and no leakage. When unfrozen, the reagents shall be clear, without sediments.
2. The consistency rates of both positive reference samples are 100%.
3. The consistency rates of both negative reference samples are 100%.
4. The kit allows the detection of 5% of specific gene mutations in a 10 ng DNA sample.
5. The repeatability is 100% by detecting designated sample for 10 consecutive times.

Precautions and Warning

1. Please read the instruction carefully in prior to experiments.
2. Conduct experiments abided by laboratory regulations to reduce cross-contaminations of products or reagents; divide experiment areas into different function zones if possible.
3. Clean experiment areas before experiment with 10% hypochlorous acid followed by water rinsing. Sterilize the environment and pipettes with 10% hypochlorous acid, 75% ethanol, or UV radiation.
4. Avoid using peripheral wells of PCR instrument; vacate holes or columns between samples to avoid cross-contamination.
5. Testing results might be influenced by sample sources, sampling process, sample quality, carriage conditions, sample handling, etc; also might it be limited by the quality of DNA, instrument types, operating environment, and the limitation of current molecular biotechnology, all of which may lead to false positive/ negative results. The users should thoroughly be informed of potential errors as well as the limitation of accuracy.
6. Avoid unnecessary freezing-thawing the reagents, the reagents were allowed to undergo no more than 5 freeze-thaw cycles.
7. The quality of DNA matters experimental results to a great extent, hence, purification of extracted DNA with magnet beads is highly suggested. Purified DNA should be stored as required ($-20\pm 5^{\circ}\text{C}$) or proceed to further steps immediately;
8. Do not substitute any original reagents contained in the kit. Do not mix reagents of different lots.
9. The use of filter tips is highly recommended to avoid false-positive results caused by contamination of reagents.
10. Be cautious of contamination from external DNA; use specific pipettes and tips for reagents preparation and template addition.
11. All reagents in use have potential hazard. For first-use of this kit, you may receive training by our technical supports. All used contents of the kit should be considered as clinical dessert and should be disposed properly.
12. All samples including positive control in the kit should be considered potential infectious substances. They should be handled carefully.

Notes

Symbol	Legend
	Indicates the need for the user to consult the instructions for use.
	Indicates a medical device that is intended to be used as an in vitro diagnostic medical device.
	Indicates the date when the medical device was manufactured.

	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Indicates the temperature limits to which the medical device can be safely exposed.
	Indicates the date after which the medical device is not to be used.
	This is the correct upright position of the distribution packages for transport or storage.
	Indicates a medical device that needs to be protected from moisture.
	Indicates a medical device that needs protection from light sources.
	Indicates the medical device manufacturer.
	Indicates the authorized representative in the European Community/European Union.
	The product meets the basic requirements of European in vitro diagnostic medical devices directive 98/79/EC.

References

1. Litton JK, Burstein HJ, Turner NC. Molecular Testing in Breast Cancer. Am Soc Clin Oncol Educ Book. 2019 Jan.
2. Wendt C, Margolin S. Identifying breast cancer susceptibility genes - a review of the genetic background in familial breast cancer. Acta Oncol. 2019 Feb;58(2):135-146.
3. Chinese Society of Pathology, Chinese Pathology Quality Control Center. Chinese expert consensus on BRCA1/2 variant interpretation. China Pathol. 2021.50(6): 565-571
4. Demir S, Tozkir H, Gurkan H, Atli EI, Yalcintepe S, Atli E, Sezer YA, Eker D, Tuncbilek N, Tastekin E, Ozen Y, Cicin I. Genetic screening results of individuals with high risk BRCA-related breast/ovarian cancer in Trakya region of Turkey. J BUON. 2020 May-Jun;25(3):1337-1347.
5. Doren A, Vecchiola A, Aguirre B, Villaseca P. Gynecological-endocrinological aspects in women carriers of BRCA1/2 gene mutations. Climacteric. 2018 Dec;21(6):529-535.
6. Tsang JYS, Tse GM. Molecular Classification of Breast Cancer. Adv Anat Pathol. 2020 Jan;27(1).
7. NCCN Clinical Practice Guidelines in Oncology: Ovarian Cancer Continue Including Fallopian Tube Cancer and Primary Peritoneal Cancer, Version 4.2022.
8. Baretta Z, Mocellin S, Goldin E, Olopade OI, Huo D. Effect of BRCA germline mutations on breast cancer prognosis: A systematic review and meta-analysis. Medicine (Baltimore). 2016 Oct;95(40).
9. Wang YA, Jian JW, Hung CF, Peng HP, Yang CF, Cheng HS, Yang AS. Germline breast cancer susceptibility gene mutations and breast cancer outcomes. BMC Cancer. 2018 Mar 22;18(1):315.
10. Edaily S, Abdel-Razeq H. Management Strategies of Breast Cancer Patients with BRCA1 and BRCA2 Pathogenic Germline Variants. Onco Targets Ther. 2022 Jul 27;15:815-826.
11. NCCN Clinical Practice Guidelines in Oncology: Breast Cancer, Version 4.2022.

Appendix Table 1. Information of 32 IIDx Recognition Sequences Based on Illumina

Strip Color	Index Number	i7 Sequence	i5 Sequence	Strip Color	Index Number	i7 Sequence	i5 Sequence
Purple	Index_001	TAAGGCGA	CTCTCTAT	White	Index_017	TAAGGCGA	GTAAGGAG
	Index_002	CGTACTAG	TATCCTCT		Index_018	CGTACTAG	ACTGCATA
	Index_003	AGGCAGAA	GTAAGGAG		Index_019	AGGCAGAA	AAGGAGTA
	Index_004	TCCTGAGC	ACTGCATA		Index_020	TCCTGAGC	CTAAGCCT
	Index_005	GGACTCCT	AAGGAGTA		Index_021	GGACTCCT	CGTCTAAT
	Index_006	TAGGCATG	CTAAGCCT		Index_022	TAGGCATG	TCTCTCCG
	Index_007	CTCTCTAC	CGTCTAAT		Index_023	CTCTCTAC	CTCTCTAT
	Index_008	CGAGGCTG	TCTCTCCG		Index_024	CGAGGCTG	TATCCTCT
Green	Index_009	TAAGGCGA	TATCCTCT	Yellow	Index_025	TAAGGCGA	ACTGCATA
	Index_010	CGTACTAG	GTAAGGAG		Index_026	CGTACTAG	AAGGAGTA
	Index_011	AGGCAGAA	ACTGCATA		Index_027	AGGCAGAA	CTAAGCCT
	Index_012	TCCTGAGC	AAGGAGTA		Index_028	TCCTGAGC	CGTCTAAT
	Index_013	GGACTCCT	CTAAGCCT		Index_029	GGACTCCT	TCTCTCCG
	Index_014	TAGGCATG	CGTCTAAT		Index_030	TAGGCATG	CTCTCTAT
	Index_015	CTCTCTAC	TCTCTCCG		Index_031	CTCTCTAC	TATCCTCT
	Index_016	CGAGGCTG	CTCTCTAT		Index_032	CGAGGCTG	GTAAGGAG



Lotus NL B.V.

Address: Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.

E-mail: peter@lotusnl.com



Manufacturer: XIAMEN SPACEGEN CO., LTD.

Address: 4th floor, No.2041 Xizhou Road, Xike Town, Tong'an District,

Xiamen 361100, P. R. China

Tel: +86 592 7578317

Fax: +86 592 7578319

E-mail: spacegen@ispacegen.com

Website: http://www.sspacegen.com